CRISTALENE / CRISTALENE MID

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INSTRUCTIONS FOR USE

IMPORTANT!

This data sheet is designed to assist you in the use of abdominal wall implants. It is not intended to be a reference document for the treatment of hernias and eventrations.

These medical devices must only be used by surgeons qualified to perform gastrointestinal surgery. In particular, the surgeon must be familiar with the surgical techniques and principles relating to the use of abdominal wall implants for the treatment of hernias and eventrations. Should you require any further information about the products, please contact the company representative.

BEFORE INSERTING THIS IMPLANT, PLEASE READ ALL OF THE FOLLOWING INFORMATION CAREFULLY.

CRISTALENE and CRISTALENE Mid are surgical implants designed to provide reinforcement in the human body. To achieve this, the prosthesis must be implanted very precisely in the tissues following a specific surgical technique. The surgeon must master this surgical technique. If the surgeon requires any further technical information or training, they must contact their supplier or ABISS-France.

WARNING

The patient must be informed of the risks / benefits of the intervention and accept the benefits and risks after a period of reflection.

These prostheses are permanent implants. The prescribing doctor must inform their patient that the use of a permanent implant may result in complications whose treatment may require a surgical procedure.

Certain complications may persist even after treatment.

PRODUCT DESCRIPTION

These prostheses are abdominal wall reinforcements available in different sizes; please contact a company representative to find out about the full range of models.

The implantable reinforcement is composed of a lock-stitch knitted mesh in polypropylene monofilament.

The mesh is knitted according to a procedure that allows it to be cut to size and to the dimensions required without unraveling.

It is a sterile, implantable, single-use device for use on one patient only.

INDICATIONS

The prostheses are indicated for reinforcement of the abdominal wall, such as cures for hernias and eventrations and other fascial defects requiring reinforcement repair.

This device can be inserted laparoscopically.

FUNCTIONING PRINCIPLE

These abdominal wall reinforcements are composed of knitted open-mesh polypropylene monofilaments, which allow tissue colonization and thus provide a solid tissue framework to maintain the prosthesis within the surrounding tissue. This ensures that the abdominal wall is reinforced and thus treats the hernia or eventration.

CONTRAINDICATIONS

All of the usual contraindications associated with the use of abdominal wall reinforcement implants also apply to these prostheses.

These include, but are not limited to:

- Pregnancy or a desire to become pregnant
- Patients who are breast-feeding
- Patients who are still growing
- Patients whose immune system is weakened or who have other disorders that could compromise healing
- Patients with an active infection or malignant tumors in the abdominal cavity
- Patients with blood clotting disorders
- Patient presenting with a hypersensitivity to synthetic derivatives such as polyamide
- Patients susceptible to alcohol or drug abuse
- Any contraindication for surgery

ADVERSE EFFECTS

The potential adverse effects are those associated with the implantation of prosthetic materials.

It is the surgeon's responsibility to inform the patient or the patient's family, prior to the operation, of the possible adverse effects associated with the use of these prostheses.

These adverse effects include, but are not limited to:

linked to the product:

- local irritation or even a foreign body reaction. Foreign body reactions may be manifest by the occurrence of an inflammatory granuloma.
- pain.
- as with all foreign bodies, the prosthesis may aggravate an existing infection.
- shrinkage leading to a defect of coverage.
- intestinal adhesions.

linked to the surgical procedure:

- the presence of seromas, hematomas, bleeding, vascular lesions, neuropathy, pain.
- the formation of a fistula, the possibility of extrusion.
- dysesthesias.
- a recurrence of the hernia or the soft tissue defect.
- sexual dysfunction can occur, including erectile dysfunction, painful ejaculation.

In all cases, the implant may need to be removed.

If the surgery is unsuccessful, it is important to fully re-evaluate the patient's case.

WARNINGS AND PRECAUTIONS FOR USE

- Before using the prosthesis, the user must be familiar with the surgical techniques and principles relating to the use of non-absorbable meshes for the treatment of hernias and eventrations.
- As talcum powder can cause inflammation, it is recommended that the prosthesis is handled with non-powdered gloves or the gloves should be rinsed thoroughly prior to handling.
- Each prosthesis must be carefully examined before the surgical procedure and undergo continuous monitoring throughout the surgical procedure to ensure that the structural integrity of the device is not compromised in any way. A damaged prosthesis must not be implanted or repaired.
- Good surgical practices must be followed when positioning the prosthesis. Standardized postoperative protocols must be followed.
- Do not implant in direct contact with the viscera or in the case of infection.
- Check the prosthesis is not folded after implantation.

- Do not let the prosthesis come into contact with potentially damaging elements.
- The suture line must not be placed less than 7 mm from the edge of the mesh.
- This device is a single-use implant. To avoid the transmission of infectious agents, do not re-use.
- After use, this product should be considered as a biological hazard. It should be handled and disposed of in accordance with accepted medical practices and current legal requirements and regulations.
- The management and postoperative monitoring of wounds and infections should be conducted in accordance with appropriate surgical practices.
- In the event of postoperative infectious complications close to the implantation site, total removal of the prosthesis is recommended.
- Before using these prostheses, it is recommended that the surgeon can expertly manage complications with synthetic implants for curing hernia and eventrations.

ADVICE TO THE PATIENT

- It has been established that perioperative smoking in the adult increases the risk of general complications (three-fold increase in the risk of infection, two to four-fold increase in the risk of healing complications). Smoking should be discontinued six to eight weeks before the procedure to lower the risk of surgical complications related to tobacco. The continuation of tobacco cessation is recommended for 3 to 4 weeks postoperatively.
- Obesity is a risk factor for infectious type complications in the postoperative period. Weight loss is advisable before any surgery requiring a reinforced repair.
- After the operation the patient should be advised to avoid significant physical strain (lifting heavy loads, sports etc.).
- The patient can resume an activity on the advice of their doctor.
- The patient should be informed to see their surgeon immediately in the event of any problem or dysfunction, particularly pain or fever.

INSTRUCTIONS FOR USE

The prosthesis is supplied in a sterile pack and must be stored at an ambient temperature.

Check that the packaging is intact prior to use. Do not use the prosthesis if the packaging is open or damaged.

Each prosthesis should be carefully examined prior to implantation. Do not use the prosthesis if it is damaged in any way.

The description below of the surgical procedure does not replace surgical training. It is the surgeon's responsibility to be trained in the surgical procedure before positioning this type of prosthesis.

The most elastic direction of the prosthesis is obtained when the two rounded corners are positioned vertically.

The surface of the implanted prosthesis must broadly cover the defect. A 4 cm overlap is recommended.

These prostheses have been developed to strengthen or close a fascial gap, as in the techniques for the surgical treatment of hernias and eventrations. Due to the different variations, we cannot describe a standard technique here. If a mesh reference is used as part of a specific surgical technique, it is essential that you refer to the literature.

The prosthesis must be held in place with fixations.

Fixation methods recommended for these prostheses are sutures, non-absorbable or absorbable staples, as well as Takers®. It is recommended that the fixation points are placed at about 7 mm from the edge of the mesh. Complications during the use of stapling systems have been described in the literature. You are advised to read the recommendations on the use of stapling systems in the case of this type of prosthesis.

TECHNICAL SPECIFICATIONS

The prosthesis is a sterile, implantable, single-use medical device for use on one patient only.

CRISTALENE: - Sizes available: 6 x 11 cm - 11 x 16 cm 7.5 x 15 cm - 15 x 15 cm 9 x 13 cm - 29 x 29 cm - Monofilaments of knitted polypropylene, non-absorbable - Diameter of the threads: 0.15 mm - Mean thickness: 0.5 mm - Mean grammage: 70 g/m² - Mean resistance to rupture > 6.56 kg/cm²	CRISTALENE MID: - Sizes available: 6 x 11 cm - 11 x 16 cm 7.5 x 15 cm - 15 x 15 cm 9 x 13 cm - 20 x 20 cm - Monofilaments of knitted polypropylene, non-absorbable - Diameter of the threads: 0.08 and 0.15 mm - Mean thickness: 0.46 mm - Mean grammage: 51 g/m² - Mean resistance to rupture > 3.08 kg/cm²
---	---

Note: Values currently measured may undergo slight variations due the limits of manufacturing procedures. Without influencing the function of the device.

EVALUATION OF THE DEVICE

ABISS asks doctors or medical staff to inform the company of any complication that follows the use of this device, and to return any explanted components or devices. To ensure safe handling when sending or receiving the devices, ABISS requests that the latter be decontaminated prior to dispatch.

SYMBOLS USED

 Manufacturer

 Do not reuse

 Sterilized using ethylene oxide

 Batch number

 Warning, read the instructions

 Do not re-sterilize

 Do not use if the packaging is damaged

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