
CYRENE

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EN

IMPORTANT!

These instructions are designed to help you use this single-use sling. This is not a reference guide for providing indications and techniques of urethral support.

Should you require any further information, please contact the company's representative.

BEFORE INSERTING THIS IMPLANT, PLEASE READ ALL OF THE FOLLOWING INFORMATION CAREFULLY.

The CYRENE sling is a surgical implant designed to provide a function in the human body. To achieve this, the CYRENE sling must be implanted precisely and meticulously in the tissues following a specific surgical technique. The physician must master this surgical technique as well as the specificities associated with the implantation of CYRENE. Standardized postoperative protocols must be followed. If the physician requires any further technical information or training, they must contact their supplier or ABISS-France.

WARNING

The patient must be informed of the risks / benefits of the operation and accept the benefits and risks after a period of reflection.

The CYRENE sling is a permanent implant. The prescribing physician must inform their patient that the use of a permanent implant may result in complications whose treatment may require a surgical procedure.

Certain complications may persist even after treatment.

PRODUCT DESCRIPTION

The CYRENE sling is composed of polypropylene prosthetic sling, whose extremities are capped with stainless steel in the premium version. This cap is used to clip the sling to the tunnelling device.

Three models are available:

- CYRENE ref. 102313: transobturator approach
- CYRENE premium ref. 112313: Out-In transobturator approach
- CYRENE premium In-Out ref. 122313: In-Out transobturator approach

INDICATIONS

The CYRENE sling is designed to provide urethral support in the surgical treatment of female stress incontinence.

CONTRAINDICATIONS

It is the surgeon's responsibility to inform the patients or their representatives of the contraindications associated with the use of this product. All of the usual contraindications that apply to the use of urethral slings also apply to the implantation of the CYRENE sling.

These include, but are not limited to, patients with the following conditions:

- Pregnancy or desire for future pregnancy
- Breast-feeding
- Potential for further growth (e.g. adolescents)
- Genitourinary infections (urinary infection, vaginal infection, etc.)
- Perioperative injuries of the urethra (fistula, diverticula)
- Pelvic radiotherapy
- Auto-immune disease, compromised immune system or other conditions that could compromise healing
- Hypersensitivity to synthetic derivatives of the polypropylene type
- Fibromyalgia

ADVERSE EVENTS

The potential adverse events are those associated with the implantation of biomaterials.

It is the surgeon's responsibility to inform the patient or her family, prior to the operation, of the possible

adverse effects associated with the use of the CYRENE sling.

These include, but are not limited to:

a) linked to the sling:

- Local irritation or even a foreign body reaction.

Foreign body reactions may take the form of inflammatory granuloma.

- Pain, some of which may become permanent, sexual dysfunction such as dyspareunia.
- As with any foreign body, the polypropylene slings can potentiate the risk of infection or can migrate in the tissue.

b) linked to the surgical procedure:

- Injury or laceration to the vessels resulting in hemorrhage, hematoma, nerve damage, or damage to the bladder, urethra, vagina, may occur during the surgical procedure, and may require surgical repair.
- Acute retention, temporary or permanent obstruction of the lower urinary tract if the sling was slung too tightly under the mid-urethra during implantation.
- Erosion, extrusion, fistula, dehiscence of the vaginal incision, failure to heal, or necrosis of the vaginal mucosa.
- Vaginal discharge or bleeding; if this occurs, the patient should consult her physician as soon as possible.
- Pain, some which may become permanent, dyspareunia.
- As with all suspension surgery, there can be de novo emergency cases, de novo cases of hyperactive bladder, and sexual dysfunction of the dyspareunia type.
- A recurrence of stress urinary incontinence.
- Neurological and/or neuromuscular symptoms.
- Vaginal extrusion or prosthetic retraction may be the cause of partner discomfort.

Certain adverse effects may disappear or resolve over time. The treatment of certain effects may require surgical removal of the CYRENE sling.

If the surgery is unsuccessful, it is important to fully reevaluate the patient's case.

WARNINGS AND PRECAUTIONS ASSOCIATED WITH USE OF THE PRODUCT

- Each CYRENE sling must be carefully examined before the surgical procedure and undergo continuous monitoring throughout the surgical procedure to ensure that the structural integrity of the device is not compromised in any way. A damaged CYRENE sling should not be implanted or repaired.
- The use of slings in patients who already have an implant or with a history of implant use, or who have undergone previous surgery in the pelvic floor area, should be evaluated with care. These patients may have a higher rate of complications.
- It is important that each procedure for the implantation of a CYRENE sling be evaluated and adapted to each individual patient. The use of imaging techniques before and after the procedure may facilitate the placement of the sling and confirm the absence of injury to adjacent anatomical structures.
- Before using this type of sling, it is recommended that the surgeon can expertly manage complica-

tions with synthetic implants for treating stress urinary incontinence.

- Good surgical practices must be followed when implanting the CYRENE sling. Standardized postoperative procedures should be followed.
- Do not let the sling come into contact with potentially damaging instruments (staples, clips or clamps).
- Check that the sling is placed without tension under the mid-urethra.
- Check the sling is not folded after implantation.
- A cystoscopy must be performed at the end of the procedure to check the integrity of the urethra and bladder.
- Only use adapted ABISS tunnelling devices.
- The caps in the premium versions are used to clip on the tunnelling devices. They cannot be used for more than 60 minutes at a time in the human body. In no event should they be implanted. **WARNING:** these specifications are only applicable for A SINGLE connection or disconnection of the tunneling device and cap. The safety of any additional use of this equipment cannot be guaranteed.
- This device is a single-use implant. To avoid the transmission of infectious agents, do not re-use.
- After use, this product should be considered as a biological hazard. It should be handled and disposed of in accordance with accepted medical practices and current legal requirements and regulations.
- The management and postoperative monitoring of wounds and infections should be conducted in accordance with appropriate surgical practices.
- In the event of post-operative infectious complications close to the implantation site, the prosthesis should be removed in its entirety.

WARNINGS AND RECOMMENDATIONS TO THE PATIENT

As with any surgical procedure, the patient must be warned about certain risk factors:

- Perioperative smoking:

It has been established that perioperative smoking in the adult increases the risk of general complications (three-fold increase in the risk of infection, two to four-fold increase in the risk of healing complications). In the event of a prosthetic surgery, the risk of vaginal erosion is markedly increased (multiplied by four). Smoking should be discontinued six to eight weeks before the procedure to lower the risk of surgical complications related to tobacco. The continuation of tobacco cessation is recommended for 3 to 4 weeks postoperatively.

- Menopause:

The surgical treatment of urinary incontinence in a menopausal patient could increase the risk of prosthetic exposure.

- Obesity:

Obesity is a risk factor for infectious type complications in the postoperative period. Weight loss is advisable before any surgery for urinary incontinence with a prosthesis.

- Diabetes
- Connective tissue disorder
- Coagulation disorder
- The patient should be informed that the beneficial effects of this surgical procedure may be lost if she becomes pregnant at any time in the future.
- The patient should be told to consult her physician immediately in the event of any problem or dysfunction, particularly in the event of dysuria, vaginal discharge, bleeding, pain or fever.
- For one month after the procedure, the patient should be warned that they must avoid physical exertion (lifting heavy loads, sport), and refrain from all sexual activity and taking baths. The patient may resume activities on the advice of her doctor.

INSTRUCTIONS FOR USE

The CYRENE sling is supplied in a sterile pack and must be stored at an ambient temperature. Check that the packaging is intact prior to use. Do not use the prosthesis if the packaging is open or damaged. The CYRENE sling can be implanted under general, regional, or local anesthesia.

The sling should be carefully examined prior to implantation. Do not use the sling if it is damaged in any way.

We will only describe the transobturator approach. Implantation via a retropubic approach is possible and can be accomplished using an appropriate tunnelling device. The description below is not the surgical procedure. It relates to points concerning the use of CYRENE. The surgeon must be trained in the surgical procedure before implanting the sling.

We will only describe the steps that are specific to the CYRENE sling.

Transobturator procedure (Out-In approach)

- After incision of the vaginal mucosa with a cold knife and opening of the fascia with Mayo scissors, on each side make 2 passages between the fascia and urethra until abutting the perineal membrane and the posterior edge of the ischiopubic rami. Do not open the perineal membrane.
- Never dissect between the fascia and mucosa as there is an increased risk of vaginal erosion, extrusion and necrosis.
- Perform the external cutaneous counter incisions of the labia majora on a horizontal line passing between the clitoris and the urethral meatus. This punctiform incision is located internally of the ischiopubic ramus.
- Passage of the instrumentation following a horizontal OUT/IN path. This passage follows 3 stages:
 - firstly, after a contact with the bone, passage of the instrumentation externally of the ischiopubic ramus
 - secondly, after the passage externally and behind the bone, sharp orientation of the instrumentation horizontally until it makes contact with the operator's finger situated between the urethra and the fascia under the ramus. Passage across the obturator muscles.

- third stage, passage of the instrumentation across the perineal membrane guided by the finger until outside the vaginal incision. Passage across the perineal membrane is performed laterally to the urethra at the level of the middle third.
- Check visually that the tunneller has not pierced the vagina.
- Introduce about 3 cm of the tapered end of the sling into the eye of the tunneller.
- Put the sling into position by reversing the tunneller. The sling must be located between the urethra and the fascia, and not between the fascia and the vaginal mucosa (risk of extrusion, erosion and necrosis).
- Ensure that the sling is not twisted under the urethra.
- It is important to leave a visible space (3 – 5 mm) between the urethra and the sling, so that the sling can be adjusted tension free.
- Closure in a single layer of the fascia and the mucosa using separate stitches of resorbable 2/0 thread with vertical mattress sutures (Blair-Donati).
- Cut off any excess sling at the level of the obturator incisions.
- Separate the sling from the skin layers, and suture with a resorbable thread suture.
- Close the vaginal incision with interrupted sutures of gradually resorbable thread.

Transobturator procedure (Out-In approach) with the premium version

- Exteriorize the tip of the tunnelling device through the vaginal incision whilst holding the finger against the tip to protect the urethra.
- Visually check that the tunnelling device has not transfixed the vagina.
- Clip the cap of the sling to the extremity of the tunnelling device, an audible click will be heard, and check that the cap is firmly attached to the tunnelling device.
- If using reusable tunnelling devices, it is advisable to use Kocher type forceps to unclip the cap and turn it through a quarter-turn to the front or back.
- Separate the sling from the cutaneous planes, and secure using one absorbable suture.

Transobturator procedure (In-Out approach) with the premium version

- Clip the cap of the sling to the helical needle, a slight click is heard, and check that it is firmly attached.
- Make a mini-incision with a scalpel to help the end of the helical needle out.
- To release the sling, either cut the sling at the level of the cap, or disconnect the cap from the helical needle (it is advisable to use Kocher-type forceps and turn it through a quarter of a turn to the front or back).
- Repeat the procedure with the second tunnelling device: clip the cap and perforate the tissues.
- To obtain correct tension in the sling, it is important to leave a visible space (3 to 5 mm) between the urethra and the sling.
- Cut any excess sling flush with the obturator incisions.
- Separate the sling from the cutaneous planes, and secure using one absorbable suture.

EVALUATION OF THE DEVICE

ABISS asks doctors or medical staff to inform the company of any complication that follows the use of this device, and to return any explanted components or devices. To ensure safe handling when sending or receiving the devices, ABISS requests that the latter be decontaminated prior to dispatch.

DESCRIPTION AND CHARACTERISTICS OF THE CYRENE SLING

The sling is a sterile, implantable, single-use medical device with the following technical characteristics:

- Monofilament woven polypropylene, non-resorbable
- Mean thickness : 0.0005 m
- Mean grammage: 0.899 g/m
- Diameter of the threads: 0.15 mm
- Size of the principal pores in mm:

* 1.31 x 0.44

* 1.62 x 1.02

- Mean porosity: 85%
- Resistance to rupture > 55 N
- Mean elongation to rupture: 58%
- Mean elongation under 10 N: 8%

Values currently measured may undergo slight variations due the limits of manufacturing procedures. Without influencing the function of the device.

The tests were performed in accordance with the norm NF S 94-801 "Vaginal reinforcement implant for stress urinary incontinence and / or pelvic organ prolapse repair surgery through vaginal approach – preclinical and clinical tests" (Dec. 2007).

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