
OPUR

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
EN

IMPORTANT!

These instructions are designed to help you use the OPUR single-use bladder hammock. This is not a reference guide for techniques of vesical support.

Before using the prosthesis, the user must be familiar with the surgical techniques and rules relative to the implantation of OPUR used for the transvaginal treatment of cystocele.

Should you require any further information about the products, please contact the company representative.

BEFORE INSERTING THIS DEVICE, PLEASE READ ALL OF THE FOLLOWING INFORMATION CAREFULLY.

OPUR is a surgical implant designed to provide a function in the human body. To achieve this, the OPUR device must be implanted very precisely in the tissues following a specific surgical technique. The surgeon must master this surgical technique as well as the specificities associated with the implantation of OPUR.

If the surgeon requires any further technical information or training, they must contact their supplier or ABISS-France.

WARNING

The patient must be informed of the risks / benefits of the intervention and accept the benefits and risks after a period of reflection.

OPUR is a permanent implant. The prescribing doctor must inform their patient that the use of a permanent implant may result in complications whose treatment may require a surgical procedure.

Certain complications may persist even after treatment.

PRODUCT DESCRIPTION

The OPUR bladder hammock kit, ref. 152000, contains 3 elements:

- One OPUR bladder hammock,
- A set of 4 nitinol wires (500 mm), used for the placement of the anterior and median arms,
- A set of 2 nitinol wires (650 mm), used for the placement of the posterior arms.

• The OPUR bladder hammock is a single-use, sterile implantable medical device, which should be used in one patient only. It is composed of woven polypropylene monofilaments.

The OPUR bladder hammock is composed of 6 arms: 2 anterior arms (blue thread), 2 median arms (green thread), and 2 posterior arms (white thread), enabling effective support.

During implantation, the blue dot on the implant should always be directed towards the top of the device.

• The nitinol wires are sterile, single-use guides. They are used to position the arms of the OPUR bladder hammock. Positioned in such a way as to form a loop in the eye of the appropriate tunneling device, the nitinol wires are used to grasp one of the straps of the bladder hammock and to guide it during placement.

• We recommend the use of OPUR ancillary instruments, ref. 942000, which are specially designed for use during the implantation of the OPUR bladder hammock. Three specific ancillary instruments (tunneling devices), corresponding to the different arms of the hammock, enable you to create the required surgical access for the placement of the hammock. The tunnellers are colour-coded with the same colour as the corresponding guide wires, which are attached to the extremities of the arms of the OPUR bladder hammock.

FUNCTIONING PRINCIPLE

The functioning principle of the bladder hammock is to reproduce the natural suspension systems of the pelvic organs according to the surgical techniques OUT/IN or IN/OUT. The anterior arms are placed in the anterior transobturator, the median arms are positioned in the posterior transobturator, and the posterior arms are placed across the sacrospinous ligament. The bladder hammock is composed of polypropylene mesh, which allows tissue colonization and thus provides a solid tissue framework to maintain the implant within the surrounding tissue.

INDICATIONS

The OPUR bladder hammock is designed to strengthen the tissues and provide long-term stabilization of the pelvic floor following genital prolapse of the anterior and middle compartments.

CONTRAINDICATIONS

It is the surgeon's responsibility to inform the patients or their representatives of the contraindications associated with the use of this product.

These include, but are not limited to:

- Pregnancy or a desire to conceive,
- Breast-feeding,
- Patients who are still growing,
- Patients with a systemic infection, infection of the urogenital tract, or infectious disorder causing altered intestinal transit,
- Immunodepressed patients,
- Patients who have had pelvic radiotherapy, major atrophy of the vaginal tissue, associated gynecological pathology,
- Patients with a clotting disorder,
- Any contraindication for surgery,
- Patient presenting with a hypersensitivity to synthetic derivatives such as polyamide.

ADVERSE EVENTS

Adverse effects following transvaginal procedures involving the insertion of prosthetic material may also occur after the implantation of the OPUR bladder hammock.

These adverse effects include, but are not limited to:

Associated with the prosthesis:

- Pain, dyspareunia.
- An inflammatory reaction, even a suppurative inflammation.
- The formation of adhesions.

Associated with the surgical procedure:

- Pain.
- The onset of erosion or the formation of a fistula, the possibility of extrusion, failure to heal.
- Vaginal discharge or bleeding; if this occurs, the patient should consult her surgeon as soon as possible.
- Injury or laceration to the vessels resulting in haemorrhage, nerve damage, or damage to the bladder, vagina, or rectum may occur during the surgical procedure, and may require surgical repair.
- Surgical repair of prolapse may be directly responsible for the onset of de novo functional disorders such as: urinary incontinence, dyspareunia, constipation, gaseous incontinence, urge incontinence, dysuria, or even decompensation of another compartment.

Certain adverse effects may disappear or resolve over time.

• As with any foreign body, the polypropylene prosthesis can increase the risk of infection.

If the procedure is unsuccessful, it is important to reassess the patient's case.

In all cases, the implant may need to be removed.

WARNINGS AND PRECAUTIONS FOR USE

- Since talcum powder can cause inflammation, the bladder hammock should be handled with non-powdered gloves, or the gloves should be rinsed thoroughly prior to handling.
- Each OPUR prosthesis must be carefully examined before the surgical procedure and undergo continuous monitoring throughout the surgical procedure to ensure that the structural integrity of the device is not compromised in any way. A damaged OPUR prosthesis should not be implanted or repaired.
- Good surgical practices should be followed when implanting the OPUR bladder hammock. Standardized postoperative procedures should be followed.
- Do not let the bladder hammock come into contact with potentially damaging instruments (clips or clamps).
- The OPUR kit comprises one implant and sterile single-use nitinol wires. To avoid the transmission of infectious agents, do not re-use.
- In no event should the nitinol wires or implantation instruments (tunneling devices ref. 942000) be used for more than 60 minutes at a time in the human body.
- The colour-coded threads facilitate the localization. In no event should they be implanted.
- After use, the elements of the OPUR kit must be treated as hazardous biological waste. They should be handled and disposed of in accordance with accepted medical practices and current guidelines and legislation.
- The management and postoperative monitoring of wounds and infections must be conducted in accordance with appropriate surgical practices.
- Given the variable morphology of each patient, it is important that each procedure for the implantation of an OPUR be evaluated and adapted to each individual patient. The use of imaging techniques before and after the procedure may facilitate the placement of the hammock and confirm the absence of injury to adjacent anatomical structures.
- In the event of postoperative infectious complications close to the implantation site, the prosthesis should be removed in its entirety.
- Before using OPUR, it is recommended that the surgeon is trained to manage complications with synthetic implants for treating prolapses.

ADVICE TO PATIENTS

- It has been established that perioperative smoking in the adult increases the risk of general complications (three-fold increase in the risk of infection, two to four-fold increase in the risk of healing complications). In the event of a prosthetic surgery, the risk of vaginal erosion is markedly increased (multiplied by four). Smoking should be discontinued six to eight weeks before the procedure to lower the risk of surgical complications related to tobacco. The continuation of tobacco cessation is recommended for 3 to 4 weeks postoperatively.
- The surgical treatment of prolapse in a menopausal patient could increase the risk of prosthetic exposure.
- Obesity is a risk factor for infectious type complications in the post-operative period. Weight loss is advisable before any surgical correction of a prolapse with a prosthesis.
- Following implantation, the patient should be advised to avoid physical strain (lifting heavy loads, partaking in sporting activities such as cycling, jogging, etc.), abstain from sexual activity, and avoid taking baths. Resumption of normal activities will be dictated by the doctor.
- After a procedure for genital prolapse, prevention of defecatory efforts is desirable.
- The patient should be told to inform her surgeon without delay in the event of any problem or dysfunction, notably if she experiences pain, bleeding, or fever.
- The patient should be informed that the beneficial effects of this surgical procedure may be lost if she becomes pregnant at any time in the future.

INSTRUCTIONS FOR USE

These instructions are recommended for the general use of this de-

vise in the treatment of prolapse of the anterior and middle compartments. The OPUR bladder hammock should only be used by doctors who are qualified for the reconstruction of the pelvic floor in general and more specifically for the OPUR device. Variations in use due to the morphology as well as the surgical history of each patient may occur.

The OPUR bladder hammock is supplied in a sterile pack and must be kept at room temperature. Check that the packaging is intact prior to use. Do not use the prosthesis if the packaging is open or damaged. Each prosthesis should be carefully examined prior to implantation. Do not use the prosthesis if it is damaged in any way.

The description below of the surgical procedure does not replace surgical training. The surgeon must be trained in the surgical procedure before implanting OPUR.

The patient is positioned in the gynaecological position with the thighs moderately flexed onto the abdomen and the buttocks slightly off the operating table. The procedure is as follows:

1 – Dissection:

- exert traction on the cervix with a pair of forceps to expose the anterior vaginal cul-de-sac.
- make a transverse vaginal incision on the anterior aspect of the cervix on the cervical side of the vaginal cul-de-sac that has been clearly exposed by traction.
- grasp the anterior vaginal section in its full thickness (vaginal wall and Halban's fascia) with 3 Allis tissue forceps, which are drawn upwards to expose the vesicovaginal plane.
- extend the vesicovaginal dissection progressively whilst everting the anterior vaginal wall.
- free up the entire inner surface of the obturator foramen and locate the ischiatic angle of the obturator foramen of the ischial spine. Release the anterior surface of the sacrosclatic ligament.

2 – Placement of the bladder hammock:

The middle of the anterior aspect of the implant is fixed using one resorbable suture to the vaginal wall in the centre opposite the trigone, above the bladder neck.

Passage of the implant arms through the eye of the tunneller and the use of the nitinol wires:

- Pass the chosen nitinol wire through the eye of the tunneller and fold it back along the length of the needle.
- Pass the tunneller through the obturator foramen. Hold the wire, which should now form a loop.
- Withdraw the tunneller. Pass the wire of the corresponding colour to the arm that is to be positioned through the loop made in the nitinol wire. Pull on the nitinol wire. Once the arm is in position, cut the end and make sure that none of the wire is left in the patient.

Passing the anterior transobturator arms (blue threads):

- Make a punctiform incision 15 mm lateral to the ischiopubic ramus horizontally with relation to the clitoral hood. Insert the tunnelling device with the blue handle into this incision.
- Slide a finger behind the obturator foramen, feel for the point of the tunnelling device and guide the device into the anterior vaginal incision.
- Using a nitinol wire and tunnelling device, pull one anterior arm of the prosthesis through the obturator foramen. Perform the same procedure on the other side.
- Position the arms without traction to prevent the prosthesis causing dysuria, as the latter is placed just above and behind the neck of the bladder.

Passing the posterior transobturator arms or middle arms (green threads):

- Insert the green-handled tunneller one finger's breadth above the ischiatic angle of the obturator foramen and on the bisection of this angle. Make a sagittal trajectory lateral to the pudendal zone marked

with the internal finger placed on the ischial spine. The needle of the tunneller comes out into the pelvic cavity lateral to and in front of the ischial spine. The extremity of the instrument is just inside the ischiopubic ramus.

- Position the arms without any traction to prevent any rope effect from the fibrotic prosthetic arm in the latero-vesicovaginal grooves.
- Fix the upper edge of the prosthesis with 2 or 3 non-resorbable sutures to the anterior aspect of the cervix.

Passing the posterior trans-sacrospinous arms (white threads):

The posterior trans-sacrospinous arms are essential for stabilising the assembly dorsally and for preventing superior recurrence of the cystocele. They are positioned using the white-handled tunneller.

- Check that the 6 nitinol wires are present along with all of the other surgical instruments used.
- Close the vaginal incision in two simple continuous planes using resorbable suture material for perfect approximation.
- Cut the excess arm; release the skin from the deep plane.

EVALUATION OF THE DEVICE

ABISS asks doctors or medical staff to inform the company of any complication that follows the use of this device, and to return any explanted components or devices. To ensure safe handling when sending or receiving the devices, ABISS requests that the latter be decontaminated prior to dispatch.

TECHNICAL SPECIFICATIONS

The OPUR prosthesis is a single-use, sterile implantable medical device with the following technical characteristics:

- Non-resorbable woven polypropylene monofilaments
- Average thickness: 0.27 mm
- Average grammage: 22 g/m²
- Wire diameter: 0.08 mm
- Average size of pores: 1.18 x 0.70 mm
- Average porosity: 59%
- Average tensile strength: 73.6 N
- Average elongation at break: 66.95%
- Average elongation under 10 N: 19.81%

Values measured at a given moment likely to show slight variations caused by the limits of the production processes. No impact on device functionality.

The tests were performed in accordance with the norm NF S 94-801 "Strengthening implants for per-vaginal placement for the treatment of stress incontinence and / or prolapse of the pelvic organs – pre-clinical and clinical trials" (Dec. 2007).

SYMBOLS USED

Manufacturer

Do not re-use

Sterile. Sterilized using ethylene oxide

Batch number

Warning, read the instructions

Do not re-sterilize


Do not use if the packaging is damaged

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